



Shanghai AnQing Medical Instrument Co., Ltd.
3&4 Floor, No.2 Building, 366 Huiqing Rd,
East Zhangjiang High-Tech Park,
201201 Shanghai, China

Tel: +86 21 61117375
Fax: +86 21 61117374
Email: global@innovexmed.com.cn



SINGLE USE FLEXIBLE BRONCHOSCOPE





**SINGLE USE STILL
NEEDS PERFORMANCE**

BRIGHT & CLEAR



HIGH IMAGE QUALITY

INNOVEX Single-Use Flexible Bronchoscope provides every patient a clean scope with high quality view - every time. Eliminate the cost and risk associated with reprocessing and cross-contamination.



Metal snake bone structure
Excellent maneuverability
Working length 620mm

Latest CMOS technology
Superior image
Field of view 110°±10%
Depth of field: 5mm-100mm



4.1mm/1.2mm



5.0mm/2.0mm



5.9mm/2.8mm

CLEAR & SHARP IMAGE

DMR-RC1-013-01 A/0