

# **SINGLE USE FLEXIBLE URETERORENOSCOPE**



SHANGHAI ANQING MEDICAL INSTRUMENT CO, LTD.  
3&4 Floor, No.2 Building, 366 Huiqing Rd,  
East Zhangjiang High-Tech Park,  
201201 Shanghai, China


Tel: +86 21 61117375  
Fax: +86 21 61117374  
Email: [global@innovexmed.com.cn](mailto:global@innovexmed.com.cn)

**SINGLE USE STILL  
NEED PERFORMANCE**

# BRIGHT & CLEAR

## HIGH IMAGE QUALITY

INNOVEX Single-Use Flexible Ureterorenoscope provides every patient a clean scope with high quality view – every time. Eliminate the cost and risk associated with reprocessing and cross-contamination.



US 31E-12-EU  
US 31E-12-US



### CLEAR&SHARP IMAGE

Dual LED illumination  
Latest CMOS technology  
Superior image  
Field of view 110°±10%  
Depth of field: 3mm-100mm



US 31D-12-EU  
US 31D-12-US